



Periodontal


Gracey Curettes	C-1
Universal Curettes	C-4
Sickle Scalers	C-6
Scalers & Curettes	C-8






Periodontal Gracey Curettes



- Gracey Curettes are used for deep scaling, root planing and periodontal debridement.
- The offset blade provides perfect working angulation for the tooth surface which means only lower cutting edge of each blade is used.



Gracey Curettes are available in the following 3 models



Standard Gracey



Note: Normally, standard models are supplied. Mention "L" for long and "S" for short in case standard model is not required.

TiN: Titanium Nitride Coating



2611-01KT
Gracey Curette Set (7 Pieces)

Contents:-

- 2611-2** Gracey 1-2
- 2611-4** Gracey 3-4
- 2611-6** Gracey 5-6
- 2611-8** Gracey 7-8
- 2611-10** Gracey 9-10
- 2611-12** Gracey 11-12
- 2611-14** Gracey 13-14
- 2611-07** 7 Instruments Cassette

Standard Gracey



2611-2 Gracey 1-2



2611-4 Gracey 3-4



2611-6 Gracey 5-6



2611-8 Gracey 7-8



2611-10 Gracey 9-10



2611-12 Gracey 11-12



2611-14 Gracey 13-14



Periodontal
Gracey Curettes



Gracey Curettes in Silver Handle



2611-02KT
Gracey Curette Set (7 Pieces)

- Contents:-**
- 2611-2SH** Gracey 1-2
 - 2611-4SH** Gracey 3-4
 - 2611-6SH** Gracey 5-6
 - 2611-8SH** Gracey 7-8
 - 2611-10SH** Gracey 9-10
 - 2611-12SH** Gracey 11-12
 - 2611-14SH** Gracey 13-14
 - 3817-07** 7 Instruments Cassette



2611-2SH Gracey 1-2



2611-4SH Gracey 3-4



2611-6SH Gracey 5-6



2611-8SH Gracey 7-8



2611-10SH Gracey 9-10




2611-12SH Gracey 11-12




2611-14SH Gracey 13-14




The blade of an universal curette has been designed with two cutting edges and round toe making it an efficient instrument to be used on all surfaces.



Langer Curettes



Columbia Curettes





Periodontal Universal Currettes



Barnhart Currettes



2612-30 BH1/2

For the side tooth area



2612-32 BH5/6

For the front tooth area



2612-32BH BH5/6

For the front tooth area

Indiana University



2612-40 IU13/14



2612-42 IU17/18



2612-42BH IU17/18

Young-Good




2612-50 GY7/8



2612-50BH GY7/8



Two cutting edges on a straight blade that end in a sharp point.



2613-10 H6/7 Anterior Premolar Hygienist



2613-12 H5/U15 Hygienist/Towner



2613-14 U15/33 Towner/Jacquette



2613-16 U15/30 Towner/Jacquette



2613-18 H5/33 Towner/Jacquette



2613-20 J30/33 Jacquette

For tooth area



2613-22 J31/32 Jacquette


For molar area



2613-24 J34/35 Jacquette



Periodontal Sickle Scalers



2613-26 J37/38 Jacquette

For Posterior

2613-28 204 Posterior

2613-30 204S Posterior

2613-32 204SD Posterior

2613-34 Siuf 204

2613-36 Spr 2/3 Spring Field

2613-38 B2/3 Ball Scaler

2613-40 T2/3 Taylor

2613-42 Ck6 Crane-Kaplan

2613-44 CI 2/3



To remove heavy calculus above or below the surface of gingival area.
The working ends come in a variety of shapes and sizes, but are always narrow.

Goldman Fox



2614-10 GF1



2614-12 GF2



2614-14 GF3



2614-16 GF4



2614-16BH GF4



2614-18 GF5 HOE



2614-20 GF6 HOE



2614-22 GF21



2614-24 GF21S/33



Periodontal Scalers & Curettes



McCall



2614-30 MC13/14S Pointed

2614-32 MC 13/14 Rounded



2614-34 MC17/18



2614-36 MTC13/14 Pointed



2614-38 MC17/18 Pointed

2614-40 MC17/18 Rounded

Cattoni



2614-50 C103/106



2614-52 C107/108



2614-52BH C107/108

Hoes



2614-60 OR6/7



2614-62 OR8/9